EN	DE	FR	ES	PT	IT
EL	NL	DA	FI	NO	SV

VOCO

Ionofil® Plus

EN **Instructions for use**
MD EU Medical device
In accordance with DIN EN ISO 9917-1

Product description:
Ionofil Plus is a radiopaque glass ionomer cement for fillings with chemical adhesion to tooth enamel and dentine. **Ionofil Plus** has a translucency similar to that of enamel.

Shades:
A1, A2, A3

Indications:
– For restoration of deciduous teeth (especially class I)
– Repair of V-shaped defects and enamel erosions as well as root erosions in the cervical area incl. class V fillings
– Cavities caused by carious lesions
– Fillings of class III and small fillings of class I (incl. extended fissure sealing) as well as build-up fillings for crowns, linings

Contraindications:
Ionofil Plus contains polyacrylic acid and parabens. **Ionofil Plus** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
Ionofil Plus is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of its intended use and the relevant product standards.

User:
Ionofil Plus should only be used by a professionally trained dental practitioner.

Use:
Preparation:
Thanks to **Ionofil Plus'** good adhesion, it is largely possible to do without additional means of retention. However, undercuts do give the filling additional hold. The layer thickness of the fillings should not be less than 1 mm. Class V cavities should be at least 1 mm deep. Avoid feather edges. In case of doubt, conventional preparation should be given preference. In case of carious defects, prepare the cavity accordingly, clean, rinse thoroughly with water and dry. No additional lining is required with a residual dentine thickness of approximately 1-1.5 mm. Pulp protection should be used in the immediate vicinity of the pulp (less than approximately 1-1.5 mm) in order to avoid irritations.

Shake powder before use.
Place the powder and liquid on a glass plate or suitable mixing pad and then mix in the powder in portions, starting small and working up to larger quantities. For 1 level measuring scoop of powder, you need 1 drop of liquid (mix ratio 4.7-5.6 g powder: 1 g liquid). Always hold the bottle vertically when dispensing the drops. A solid plastic spatula should be used for the mixing.

Mixing time	30-40 s
Working time (15-23°C)	2-3 min.
Setting time (37°C)	2.5 min.
Total setting time	5-6 min.

The mixing time is 30-40 s (it should not exceed 1 min). The working time **following** the mixing is approximately 2-3 min. (the temperature should be 15-23°C). After its introduction into the cavity, **Ionofil Plus** sets within approximately 2.5 min.

Application:
Keep the operating site dry when placing the filling. Glass ionomer cements are sensitive to moisture (breath, saliva) during the setting phase and particularly directly after their application. For this reason, if possible, apply a protective varnish immediately. Once **Ionofil Plus** has set, 5-6 min. after application at the earliest (depending on mixing temperature and ratio), any excess material can be removed and the finishing performed (flexible discs). Then coat with varnish again. The best finishing results are obtained approximately 24 hours after placement of the filling.

Warnings, precautionary measures:
– Reclose the powder and liquid carefully after use so as to prevent moisture from entering and preserve the excellent strength. Avoid contact with mucous membranes and eyes.
– Please also do not use any insulating material for the dentine in order to preserve full adhesion (with the exception of pulp protection).
– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
Fluorosilicate glass, polyacrylic acid, water, tartaric acid, glucono delta-lactone, parabens, colour pigments

Storage instructions and application method:
Storage at 4-23°C. Reclose powder and liquid containers well immediately after use. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Ionofil Plus** must be reported to VOCO GmbH and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of **Ionofil Plus** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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