

NEEDLE HOLDER

Used to grasp and guide the needle when suturing



Mayo-Hegar Straight 16.0cm
NHMH



Mayo-Hegar Curved 16.0cm
NHMHC



Crile-Wood 15.0cm
NHCW



Adson Straight 17.0cm
NHA



Needle Holder 18.0cm
NHMH18
(Also available in 20.0cm #NHMH20)



Olsen-Hegar Needle Holder
with Scissors 17.0cm
NHOH



Derf Straight 12.5cm
NHD



Derf Curved 12.5cm
NHDC

Castroviejo
Needle Holder



14.0cm
Curved
NHCVC



14.0cm
Straight
NHCVS

Backhaus Towel Clamp



Backhaus 9.0cm
TC3



Backhaus 11.0cm
TC4



Backhaus 13.0cm
TC5

NEEDLE HOLDER TUNGSTEN CARBIDE

Needle holder feature tungsten carbide insert with a ground pyramidal profile to assure firm, positive needle grip. Rotation of needle and /or slipping of suture material is virtually eliminated



Mayo-Hegar 16.0cm
NH5042



Crile-Wood 15.0cm
NH5038



Halsey 13.0cm
NH5037



Derf 11.5cm
NH5032

Micro Castroviejo Needle Holder Tungsten Carbide Insert



Castroviejo TC Straight 14.0cm
NH5020



Castroviejo TC Curved 14.0cm
NH5021



Castroviejo TC Straight 18.0cm
NH5024R

Sterilising Forceps



Cheate Forcep 27.0cm
CF

Dressing Forceps



Gross Maier Straight 20.0cm
SDFS



Gross Maier Curved 20.0cm
SDFC



Foerster Straight 20.0cm
SDFF



Rogge 20.0cm / 30.0cm
RF20
RF30