

Prime&Bond universal™ Universal Adhesive



Product Description

Prime&Bond universal provides well balanced hydrophobic and hydrophilic properties to ensure complete coverage and penetration at varying moisture levels. So even if dentin is overly wet or dry, you get a reliable bond that actively guards against the leading causes of adhesive failure. Prime&Bond universal provides:

- Versatility for use with all etching methods, all indications
- Reliable performance on over-wet and over-dried dentin
- Virtually no post-operative sensitivity
- Low film thickness for a reduced risk of pooling and misfit of the restoration
- A relaxed working time and reduced waste, remaining usable for 30 minutes in a closed CliXdish
- High initial bond strength with Calibra Ceram cement for indirect restorations
- A patented formulation free of HEMA, TEGDMA, or bisphenol

Indications

- Direct, light-cured composite and compomer restorations.
- Composite, ceramic and amalgam repairs.
- Cavity varnish for use with fresh amalgam.
- Indirect restorations cemented with Calibra® Universal and Calibra® Ceram.

In combination with Dentsply Sirona Self Cure Activator:

- Indirect restorations cemented with dual-cure/self-cure resin cements.
- Direct dual-cure/self-cure composite restorations and core build up.



Flip-top cap for one-handed opening



Suitable for all etching techniques: Self-Etch, Selective Enamel Etching, Etch&Rinse



Prime&Bond universal remains usable up to 30 minutes in a closed CliXdish.

Features

Benefits

Active moisture control	Reliable bonding performance, no matter if the dentin is still too wet or already too dry - for virtually no post-operative sensitivity.
Active cavity spreading	Actively spreads in the cavity and helps to cover all cavity surfaces.
Self levelling	The adhesives self-levels to form a smooth and uniform adhesive layer, helping to avoid dry spots or adhesive pooling.
Low film thickness	Adhesive pools can show up as translucent areas on a radiograph, which can easily be misdiagnosed as a void. Thick adhesives may also change crown prep geometry impacting crown fit.

Prime&Bond universal™ Universal Adhesive

Prime&Bond universal Bottle

Prime&Bond universal Bottle Standard Refill

1 x 4 ml Bottle Prime&Bond universal
1 x Instructions for Use
1 x Illustrated Technique Guide

60667350



Prime&Bond universal Bottle Mini Refill

1 x 2.5 ml Bottle Prime&Bond universal
1 x Instructions for Use
1 x Illustrated Technique Guide

60667351



Self Cure Activator

For use in conjunction with Prime&Bond universal, Prime&Bond XP, Prime&Bond NT, prime&bond one ETCH&RINSE or Spectrum bond

Product Description

Self Cure Activator (SCA) is for use with Prime&Bond XP, Prime&Bond NT, Prime&Bond universal, prime&bond one ETCH&RINSE or Spectrum bond in conjunction with a dual-cure or self-cure resin cement or restorative such as Calibra, SmartCem2 or core-X flow.

SCA combined with the adhesives mentioned above is suitable for bonding of all indirect restorations including metal, porcelain and composite based crowns, inlays, onlays, veneers, bridge retainers and endodontic posts.

SCA Bottle

SCA Bottle Refill

1 x 4.5 ml Self Cure Activator
1 x Instructions for Use

634354K

