

**BURNISHER**

Designed to condense, smooth, carve and polish amalgam.



**BB21B**  
#1 #3



**BB21**  
#1 #3



**BBL3**  
#1 #3



**PKT3**  
#1 #3  
#4



**BB22/**  
23  
#1 #3 #4



**BB18**  
#1 #3



**BB27/**  
29  
#1 #3




**TB2**  
#1 #3

**PLACEMENT INSTRUMENTS**

Used to place base or liner within cavity preparations.



**PICH8**  
#1



**PICH6**  
#1 #3

Calcium hydroxide or glass ionomer base/liner placement instrument. Also useful as a small burnisher.



**SP60616**  
#1 #3

Calcium hydroxide or glass ionomer base/liner placement instrument combined with a mini-spatula for efficient mixing.



**PICH7**  
#1 #3

Calcium hydroxide or glass ionomer base/liner placement instrument.