

Bite Registration Wax



Thickness: 1.7mm to 2.0mm (Custom Thickness Catered)

Size: 160mm X 90 mm

Packing: Each Packet of 'MAARC' Bite Registration Wax contains 10 sheets wrapped in non-adhesive Silicone Paper to facilitate easy removal. Available in 250 g pack approx.

Bite Wafers pack of 10

Available in: 3 variants - Hard, Extra Hard & Zero Tolerance

Packing: Bite Registration Wax (Aluminium Filled)

Product Code: 2202/010

Packing: Classic BRW Extra Hard (Aluminium Filled)

Product Code: 2207/010

Packing: Perfect Bite (Zero tolerance)

Product Code: 2213/010

Packing: Bite Wafers

Product Code: 2210/010

Properties & Advantages:

- Contains Aluminium particles which gives even heating
- Dimensionally stable, does not shrink or distort.
- Ideal for checking and maintaining accurate occlusal records.
- When heated at 55°C it turns in to 'Dead Soft' in less than 1 minute.
- Adapts to every details and sets by water or room temperature i.e. 27°C

Benefits:

- **Accurate Bite:** Even Heating Due to Aluminum Particles.
- **No Shrinkage:** Negligible Shrinkage Due to Crystallization around Aluminum & Resin.
- **No Distortion:** Very Tough so Lesser Change of Distortion

Shelf life and storage: Best before 60 months under proper storage.