

DENTAL SYRINGE

GDC syringes are available in two styles, the cooke-waite and with large diameter finger rings for enhanced control, large window for enhanced carpule visibility and sharp harpoons for easy cartridge engagement.

Dental Syringe with colour coded ring



SYRACC1.8R
(Also available in different colours of ring)



Self Aspirating Syringe



SYRSA1.8

Intraligamental Syringe - Inject



SYRI1

"Suitable for upper and lowers, just turning the front barrel"



Not Aspirating

SYR1.8
SYR2.2

Aspirating

SYRA1.8
SYRA2.2

Aspirating

SYRF1.8

(Also available in European & American thread needle type)

ADVANTAGES OF INTRALIGAMENTAL ANESTHESIA :

- Automatic limitation of the power pressure.
- Anesthetizing of a single tooth.
- No numbness of the tongue, lips or cheek.
- Simultaneous treatment of teeth in different quadrants.
- Supports regional block anesthesia.
- Eliminates the risk of haematoma.
- Lower emotional impact on patient.

WHY INTRALIGAMENTAL IS BETTER

Because it's chrome free.
Because you can get a better sterilization.
Because you can use for long time
Because has anti-reflection surface.

OPERATING FEATURES:

- Fully autoclavable.
- Removable headpiece with angled threaded nozzle for the injection needle.
- Viewing chamber with safety sleeve for carpule cartridge.
- Handle servo dosing lever with optimal support for index finger for a gradual advance of the dosing plunger resetting key.
- No waiting time between injections and beginning of treatment.
- Excellent depth of anesthesia combined with shorter duration.
- Much less pain than standard injections (infiltration anesthesia in hard palate, intrapulpal or intraseptal injections or for papillectomy).
- A smaller quantity of given anaesthetic (~0,06 ml for single injection/click.)
- Needle suggested: 30G.