

EXTRACTION FORCEPS
AMERICAN PATTERN




American Pattern

Suggested Pair




53R# Molars
First & Second Upper
Molar Teeth Right

F53R



53L# Molars
First & Second Upper
Molar Teeth Left

F53L




Suggested Pair




88R# Nevius
First & Second Upper
Molar Teeth Right

F88R




88L# Nevius
First & Second Upper
Molar Teeth Left

F88L



23# Cowhorn
First & Second
Lower Molars

F23




Peedo Universal




150# Cryer Universal
Upper Incisors, Canines,
Premolars & Roots

F150




151# Cryer Universal
Lower Incisors, Canines,
Premolars & Roots

F151




69# Tomes
Upper & Lower
Fragments & Roots

F69




150S# Peedo Universal
Upper Primary Teeth
& Roots

F150S



151S# Peedo Universal
Lower Primary Teeth
& Roots

F151S



DID YOU KNOW ?

All instruments should be cleaned and sterilised in an open position.