

Plus  
**Metapaste**  
Water-based Calcium Hydroxide Paste

Premixed Water-based Calcium Hydroxide Paste

# Metapaste Plus



2g Syringe  
Disposable tips

Plus  
**Metapaste**  
Water-based Calcium Hydroxide Paste

### Faster Dissociation

Faster medicinal effect & easily removed

### PPG Free

Reducing Post-op sensitivity, ensuring more stable treatment

### Antibacterial with high pH: 12.15 (1day) / 12.35 (14day)

Increasing success probability of root canal treatment due to the high antibacterial effect

### Improved flowability: 19-20mm (ISO 6876:2012)

Convenient and accurate placement to the apex enhanced

### Premixed Syringe type

Quick and Easy treatment with minimal waste & clean up

## Advantages

### Applications

- Weeping Canal
- Formation of a hard tissue barrier (Apexification)
- Temporary root canal filling
- Internal & external root resorption

### Composition

- Calcium Hydroxide, Barium Sulfate, Purified Water

### Shelf life

- Two years from the production date

### Components

- 2g syringe 1EA, Disposable Tips 20EA
- 2g syringe 2EA, Disposable Tips 20EA

Plus  
**Metapex**  
Calcium Hydroxide Temporary Filling Material



Premixed Calcium Hydroxide Paste with Iodoform

# Metapex Plus

## Antibacterial with high pH: 12.00 (1day) / 12.4 (14day)

Increasing success probability of root canal treatment due to the high antibacterial effect

## Improved flowability: 19-21mm (ISO 6876:2012)

Convenient and accurate placement to the apex enhanced

## Premixed syringe type

Quick and Easy treatment with minimal waste & clean up

## High radiopacity: 5.4 Almm

Accurate diagnosis available

## Advantages

### Applications

- Weeping Canal
- Formation of a hard tissue barrier (Apexification)
- Root canal filling for deciduous teeth
- Internal & external root resorption

### Composition

- Calcium Hydroxide, Iodoform, Silicon oil

### Shelf life

- Two years from the production date

### Components

- 2.2g syringe 1EA, Disposable Tips 20EA
- 2.2g syringe 2EA, Disposable Tips 20EA

#### Meta Biomed Co., Ltd - Head Office

270 Osongsaengmyeong1-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungbuk, Korea

Phone: +82-43 216 0433  
dental@meta-biomed.com  
www.meta-biomed.com

#### Meta Biomed Inc.

3015 Advance Lane Colmar, PA 18915, USA

Phone: +1 267 282 5893  
america@metabiomed-inc.com  
www.metabiomedamericas.com

#### Meta Biomed Europe GmbH

Wiesenstr. 35, 45473 Mülheim an der Ruhr, Germany

Phone: +49 208 309 91910  
europe@metabiomed-inc.com  
www.meta-europe.com